alspl

Changing Lives For The Better....

NUCATH™

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Instructions For Use**Safety Considerations:**

Do not use if package is open or damaged in any way.

Please read all packaging insert warnings, precautions, and instructions before use. Failure to do so may result in severe patient injury or death.

This product is only designed, intended for use as indicated.

Please do not alter the catheter or any component in any way during the procedure.

Rx Only. The catheter must be used by trained and experienced personnel, who are well versed in catheterization techniques, anatomical landmarks, and potential complications.

The product is supplied sterile by Ethylene Oxide gas. Do not reuse to re-sterilize the catheter in any event.

The product must be stored in a dry environment, away from sunlight.

Description:

The NUCATH™ is a wedge pressure catheter with a balloon at its distal tip and is intended to be used in right heart catheterisations for the measurement of pressures. The lumen with the longer blue hub is meant for balloon inflation / deflation, and the shorter black hub is for the catheter to track over a guidewire.

The catheter is supplied as a kit, consisting of a 3 ml Luer-Lock syringe, one-way stopcock, and the catheter itself. It is highly recommended to only use the supplied syringe for balloon inflation.

The catheter has markings at every 10 cm interval from the distal tip. There is a ± 1 cm tolerance on the first mark. The balloon inflates to the stated diameter at the maximum inflation volume as indicated on Table 1. The diameter of the balloon is $\pm 10\%$ at the maximum inflation volume. The intended patient population comprises paediatrics and adults. In individual cases, the treating cardiologist determines the appropriate patient population in consideration of the contraindications, the health condition of the patient and based on the experience gained during the training as a cardiologist.

Indications for Use:

The NUCATH™ Wedge Pressure Catheter is indicated for measuring pressure in the right heart (including central venous pressure, right ventricle pressure, pulmonary artery pressure, and pulmonary artery wedge pressure).

Contraindications:

There are no specific contraindications for the wedge pressure catheter, other than the standard risks associated with the insertion of a cardiovascular catheter.

Warnings and Precautions:

1. Warning: Sterile, Single Use: Do not reuse, reprocess or resterilise under any circumstance. Doing so creates a potential risk of serious injury, infection, or balloon rupture, which may lead to death of the patient.
2. Warning: The balloon MUST always be inflated using Carbon Dioxide.
3. Warning: Do not exceed the maximum inflation volume mentioned on the product label. Over inflation will greatly increase the possibility of balloon rupture.
4. Warning: The balloon must always be deflated fully before insertion and removal of the catheter into the introducer sheath.

5. Precaution: Do NOT use liquids as balloon inflation media under any circumstance. This reduces the ability of the catheter to flow with blood circulation.
6. Precaution: It is highly recommended to only use the syringe provided with the catheter to inflate the balloon. In any event when the syringe supplied with the catheter is unavailable, use a luer lock 3ml syringe.
7. Warning: Transient obstruction of the pulmonary outflow tract, balloon rupture without serious sequelae, premature ventricular contractions, inter-luminal leaks resulting in air embolism and knotting of the catheter have been reported as complications of the use of balloon catheters in cardiac diagnosis and monitoring.
8. Warning: Cardiac perforation is a known complication of paediatric cardiac catheterization.
9. Warning: Care must be exercised when passing the catheter in patients with left bundle branch block because right branch block induced by traumatic catheter passage could result in complete heart block and asystole.
10. Precaution: The catheter must be carefully and gently used with fluoroscopic guidance to facilitate the complex movement within the heart.
11. Precaution: An introducer sheath and guidewire, of the recommended size (mentioned on product label) should be used with the catheter. A smaller introducer sheath, than the recommended size, must never be used as it could potentially damage the balloon. A thicker guidewire than the recommended size must not be used, as it will get stuck in the catheter.
12. In the event that there is a serious incident that has occurred in relation to the device, please report such an occurrence to the company (contact details at the bottom), and the competent authority in your state.

Description of the Procedure:

Testing the balloon before use:

1. Remove the balloon gently from the packaging. Make sure the balloon is not damaged in any way.
2. Open the stopcock, pull a vacuum with a luer lock syringe, and close the stopcock.
3. Fill the syringe with the recommended volume of inflation media for the size of catheter being used. Do not, under any circumstance, inflate more than the maximum volume as stated on the product label.
4. Attach the syringe to the stopcock.
5. Open the stopcock, and inject the inflation media, and close the stopcock. Precaution: If the balloon deflates is less than 1 minute, a leak is present in the catheter, or the stopcock, and the catheter MUST NOT be used.
6. Once it is confirmed that there is no leak present, deflate the balloon, by opening the stopcock.

Catheter insertion:

7. Flush the guidewire lumen (black hub) of the catheter. If a guidewire is already in place inside the patient, insert the catheter over the guidewire from its distal end. If there is no guidewire in place, insert an appropriate size guidewire into the guidewire lumen (black hub) from the proximal end until it comes out from the distal end of the catheter.
8. Open the stopcock, pull a vacuum on the balloon with a luer lock syringe, and close the stopcock.
9. Insert the catheter through an appropriate size introducer sheath.
10. Advance the catheter under fluoroscopic guidance until the tip of the catheter is in the right atrium.
11. Inflate the balloon so the catheter tip can be carried by the blood flow into the chambers and vessels as desired. The catheter must be manipulated very carefully to avoid kinking or knotting. Precaution: Do not introduce more than 10 cm of the catheter into the right atrium or right ventricle if the tip of the catheter is still in that chamber. To avoid knotting in such an event, deflate the balloon, withdraw the catheter from the chamber, re-inflate the balloon, and advance again. Use the markings on the catheter shaft for guidance if needed.

12. Upon completion of the procedure, deflate the balloon and remove the catheter.

Disposal after Use:

After use, medical products and accessories pose a potential biological hazard. For this reason, the products and their accessories should be handled and disposed in accordance with recognised medical procedure and in compliance with the relevant legal regulations and local ordinances.

Selection of the NUCATH™ Wedge Pressure Catheter:

The following table contains information allowing the healthcare professional to verify if the device is suitable and to select the corresponding accessories:

Size	Guidewire	Introducer Sheath	Maximum Inflation Volume	Balloon Diameter	Reference Number
4 Fr	0.021’’	5 Fr	1.5 ml	6.5 mm	NWP 4F
5 Fr	0.025’’	6 Fr	1.5 ml	8 mm	NWP 5F
6 Fr	0.035’’	7 Fr	2.5 ml	10 mm	NWP 6F

Potential Complications / Adverse Effects:

Potential balloon separation following balloon rupture or misuse may lead to the subsequent need to use a snare or other medical interventional techniques to retrieve the pieces.

Potential complications and risks relating to right heart catheterisation include:

- Hematoma at puncture site which may cause bruising around the area where the catheter was inserted
- Injury to the vein
- Cardiac perforation

Very rare complications include:

- Cardiac arrhythmias (including tachycardia) requiring treatment
- Cardiac tamponade
- Embolism caused by blood clots at the tip of the catheter
- Infection
- Low blood pressure
- Transient obstruction of the pulmonary outflow tract
- Balloon rupture without serious sequelae
- Premature ventricular contractions
- Inter-luminal leaks resulting in air embolism
- Knotting of the catheter have been reported
- Heart block
- Hypoxia
- Pulmonary haemorrhage

Warranty and Limitations:

Catheters and accessories are sold in ‘as is’ condition. The entire risk as to the quality and performance of the catheter is with the buyer. ALSPL disclaims all warranties, expressed or implied, with respect to catheters and accessories, including but not limited to, any implied warranty of merchantability or fitness for a particular purpose. ALSPL shall not be liable to any person for any medical expenses or any direct or consequential damages resulting from the use of any catheter or accessory or caused by any defect, failure, or malfunction of any catheter or accessory, whether a claim for such damages is based upon warranty, contract, tort, or otherwise. No person has any authority to bind ALSPL to any representation or warranty with respect to catheters and accessories.

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References:

1. Bierman, Howard R. "On the Intravascular Knotting of Catheters." *Vascular Surgery*, vol. 6, no. 4, 1 Sept. 1972, pp. 155–158., doi:10.1177/153857447200600401
2. Burgener, F A, et al. "Complications and Hazards with Angiographic Occlusion Balloon Catheters." *Radiology*, vol. 140, no. 3, Sept. 1981, pp. 647–650., doi:10.1148/radiology.140.3.7280230.
3. Civetta, J M, and J C Gabel. "Flow-Directed Pulmonary Artery Catheterization in Surgical Patients: Indications and Modifications of Technics." *Annals of Surgery*, vol. 176, no. 6, Dec. 1972, pp. 753–756.
4. Kelly, David T., et al. "Double-Lumen Flotation Catheter for Use in Complex Congenital Cardiac Anomalies." *American Heart Association*, vol. 44, Nov. 1971, pp. 910–913., doi:10.1161/01.CIR.44.5.910.
5. Stanger, Paul, et al. "Use of the Swan-Ganz Catheter in Cardiac Catheterization of Infants and Children." *American Heart Journal*, vol. 83, no. 6, June 1972, pp. 749–754., doi:10.1016/0002-8703(72)90205-0.
6. Swan, H. J. C., et al. "Catheterization of the Heart in Man with Use of a Flow-Directed Balloon-Tipped Catheter." *New England Journal of Medicine*, vol. 283, no. 9, 1970, pp. 447–451.
7. "Swan-Ganz - Right Heart Catheterization." Edited by Michael A Chen and David Zieve, *MedlinePlus Medical Encyclopedia*, U.S. National Library of Medicine, 2 Aug. 2016, medlineplus.gov/ency/article/003870.htm.
8. ISO 15223-1:2023, Medical Devices- Symbols to be used with medical device labels, labelling, and information to be supplied- part 1: General requirements.

Manufactured By:
Advanced LifeSciences Pvt. Ltd.
D-22 Okhla Industrial Area Phase – 1,
New Delhi – 110020, India.
Website : www.alspl.com,
E-mail: customerservice@alspl.com
Phone: +91 9818237529

IFU/07/01

Revision: 04

24th July 2024