
alspl

AD-BAND™

Radial Artery Compression Device
Dispositif de Compression de l'Artère
Radiale
Dispositivo de Compressão da Artéria
Radial
Dispositivo de Compresión de la Arteria
Radial

Table of Contents

ENGLISH (en)	2
FRENCH (fr)	3
PORTUGUESE (pt)	5
SPANISH (es)	7

Instructions for Use

Description:

Radial Artery Compression Device is principally used for compression of the radial artery after an interventional procedure to create hemostasis. The device is a band with an inflatable balloon which compresses the artery.

The device is made of the medical high molecular materials that are commonly used globally such as polyurethane. This product can apply mechanical pressure at the vessel puncture site in place of hand compression hemostasis so as to make the patients more comfortable and reduce the workload of medical staff.

Indications for Use:

ALSPL's Radial Artery Compression device is a device used for compression to create haemostasis after radial interventional procedures.

Contraindications:

- Patients hypersensitive to the materials of compression device.
- Patients whose wrist diameter at the site of puncture is too large or too small, which exceeds the allowable range of compression device.
- Patients with infection or other serious skin diseases at the site of puncture.
- As Radial Artery Compression Devices are only an accessory for procedures, the contraindications are limited to the contraindications for the respective procedures

Warnings and Precautions:

1. Please read all packaging insert warnings, precautions, and instructions before use. Failure to do so may result in severe patient injury or death.
2. Warning: Sterile by EO, Single Use: Do not reuse, reprocess or re-sterilize under any circumstance. Blood and other foreign material can collect in the inner parts of the product which cannot be removed.
3. Warning: Do not use if the packaging is damaged or opened unintentionally before use. This will lead to loss of sterility.
4. Warning: Rx Only. The products must be used by trained physicians.
5. Precaution: Aseptic techniques must be used wherever possible, especially during removal from packaging.
6. Precaution: Do not put excessive force or load on the balloon as it could lead to a leak or burst.
7. Precaution: Do not overinflate the band as it may leak or burst with over inflation.
8. Precaution: If there is itching or redness of the skin while compression, stop using the device and treat appropriately.
9. Precaution: Always check for hemostasis before removal of the device.

Instructions for Use:

1. Withdraw the introducer sheath for 2-3 cm after procedure.
2. Align the centre site of the compression device with the puncture site, and fix the fastener onto the wrist. For the tightness, it is best to fix the compression device up on the wrist so it can barely move.
3. Inject an appropriate volume of air into the balloon through the unilateral valve with the inflation device in the package in order to inflate the balloon.
4. Single balloon type: The nominal injection volume of air for is 15mL and the maximal injection volume of air is 18ml.
5. Remove the sheath, and observe whether there is haemorrhage at the site of puncture. If there is haemorrhage, inject more air (the total volume should not exceed 18mL until haemorrhage stops.)
6. Examine the haemostatic conditions, and adjust the air pressure with the inflation device.

Note: If the fastener on the compression is released before permanent haemostasis is achieved at the radial artery puncture site, haemorrhage may occur. Accidental release of the compression device could be the result of inappropriate fastener application, inappropriate compression device size selection or inappropriate compression device manipulation.

7. Before removing the compression device, confirm that appropriate time has passed in accordance with doctor's instructions and permanent haemostasis has been achieved.

Note: The injection volume of air and compression time may vary with the conditions of patients, usage and the dose of the anticoagulants i.e. heparin and the size of puncture site. Keep on observation on the puncture site and adjust the compression pressure accordingly. Avoid long-time high- pressure compression (e.g. greater than 12 hours or per doctor instructions/protocol).

Disposal:

1. Radial Artery Compression Devices can be a biological hazard if disposed incorrectly. They should be disposed using recognised medical procedures and industry best practices. Local laws and regulations should be taken into account during disposal

Potential Complications / Adverse Effects

Itching, or redness of skin. Radial Artery Compression Devices are accessories used in various medical procedures. The potential complications or adverse effects caused during use are limited to the procedure being performed.

Warranty and Limitations

Radial Artery Compression Device are sold in 'as is' condition. The entire risk as to the quality and performance of the product is with the buyer. ALSPL disclaims all warranties, expressed or implied, with respect to its products, including but not limited to, any implied warranty of merchantability or fitness for a particular purpose. ALSPL shall not be liable to any person for any medical expenses or any direct or consequential damages resulting from the use of any product or caused by any defect, failure, or malfunction of any product, whether a claim for such damages is based upon warranty, contract, tort, or otherwise. No person has any authority to bind ALSPL to any representation or warranty with respect to its products.

Notice: Any serious incident which occurs in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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Advanced LifeSciences Pvt. Ltd.