**EASY™, INDEFLATE™, DIVINE™**

**INFLATION DEVICE  
Dispositifs de Gonflage  
Aparelhos de Inflação  
Dispositivos de Inflación**

**Table of Contents**

ENGLISH (en) .....	2
FRENCH (fr) .....	3
PORTUGUESE (pt) .....	5
SPANISH (es) .....	7

## Instructions for Use:

### Description:

ALSPL Inflation Devices are disposable intervention accessories which allow the user to inflate balloon catheters during PTA, PTCA or any other procedures. The devices have a pressure gauge on for the ease of the user. There is a tubing attached to a 3-way stopcock with male luer lock for connection to catheters. The inflation devices have different types as preferred by the user.

### Indications for Use:

ALSPL's Inflation Devices are intended to be used while performing balloon dilation. The devices can be used to inflate, monitor the pressure and deflate the balloon.

### Contraindications:

- As Inflation Devices are only an accessory for procedures, the contraindications are limited to the contraindications for the respective procedures

### Warnings and Precautions:

1. Please read all packaging insert warnings, precautions, and instructions before use. Failure to do so may result in loss of blood and delay to the procedure.
2. Warning: Sterile by EO, Single Use: Do not reuse, reprocess or re-sterilize under any circumstance. Blood and other foreign material can collect in the inner parts of the product which cannot be removed.
3. Warning: Do not use if the packaging is damaged or opened unintentionally before use. This will lead to loss of sterility.
4. Warning: Rx Only. The products must be used by trained physicians
5. Warning: Only liquid inflation media such as contrast media must be used as recommended by the balloon catheter manufacturer
6. Warning: Maximum pressure should not exceed the balloon catheter manufacturer's recommendation, or 30 atm, whichever is lower
7. Precaution: Aseptic techniques must be used wherever possible, especially during removal from packaging.
8. Precaution: Connections should only be made with devices compliant to ISO 80369-7
9. Precaution: Ensure connector tubing is completely free of air before use

### Instructions for Use:

#### Preparation of the Inflation Device

1. Open the Inflation Device aseptically from the packaging
2. All aspiration (pull) or injection (push) actions should be done when the device piston is in the unlock position, lock lever on the left side for DIVINE and INDEFLATE models. For the EASY model, the trigger must be pressed while pulling or pushing to unlock the piston
3. Fill the syringe with the contrast media. Ensure there is no air in the system by purging. The device can be tapped or shaken lightly to get air to the top. If a 3-way stopcock is used, ensure it too has no air trapped inside

#### Attachment of the Inflation Device to the Balloon Catheter:

4. Prepare the balloon catheter as required by its manufacturer
5. Lock the piston of the device. This is done by turning the lever to the right in DIVINE and INDEFLATE models and releasing the trigger in the EASY model
6. Create a fluid-fluid connection between the inflation device / 3-way stopcock and catheter and hand tighten to create a leak free connection

## Use of Inflation Device

7. Rotate the knob clockwise to inflate or increase pressure and anticlockwise to deflate or decrease pressure
8. Once the dilation is complete deflate balloon but rotating the knob, or unlocking and pulling the piston

### Disposal:

9. Inflation Devices can be a biological hazard if disposed incorrectly. They should be disposed using recognised medical procedures and industry best practices. Local laws and regulations should be taken into account during disposal

## Potential Complications / Adverse Effects

Inflation Devices are accessories used in various medical procedures. The potential complications or adverse effects caused during use are limited to the procedure being performed.

## Warranty and Limitations

Inflation Devices are sold in 'as is' condition. The entire risk as to the quality and performance of the product is with the buyer. ALSPL disclaims all warranties, expressed or implied, with respect to its products, including but not limited to, any implied warranty of merchantability or fitness for a particular purpose. ALSPL shall not be liable to any person for any medical expenses or any direct or consequential damages resulting from the use of any product or caused by any defect, failure, or malfunction of any product, whether a claim for such damages is based upon warranty, contract, tort, or otherwise. No person has any authority to bind ALSPL to any representation or warranty with respect to its products.

**Notice:** Any serious incident which occurs in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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