 <small>Changing Lives For The Better...</small>	<b>Advanced LifeSciences Pvt. Ltd.</b>	<b>Document No.-</b> QSP/SOP/051		
		<b>Page No.</b>	Page 1 of 22	
	<b>Product Information Release</b>		<b>Effective Date</b>	01/06/2023
			<b>Revision No.</b>	01
			<b>Applicable Document:</b>	QSP/SOP/05

The form should be used to share the following product information. Please tick as appropriate:

**Labeling Information:**

Information that is applied to the product and to the packaging or attached to the product.  
(eg product imprint, packaging layout, label, instructions for use)

**Marketing Materials:**

Product and application information for the customer / user  
(eg in the form of brochures, brochures, flyers, in printed or electronic form (eg content for homepage))

**Patient Information:**

Special product and application information intended for the patient  
(eg patient guide for implants, patient passport)

**Technical Information:**

Special product and application information intended for the physician / user.  
(eg instructions for use / procedure)

**Training Information:**

Product information for the training of medical device consultants and users at ALSPL

Description	Instruction for use (IFU) for Guidewires
-------------	--

Notes on printing: (in printed matter)	A4 Booklet Printing.
---	----------------------


Document No: IFU/14/01	Revision No:05
------------------------	----------------

**Document History:**

Sr. No.	Date	Revision No.	Description of Revision/Change
1.	20/03/2020	00	New Document Created
2.	02/09/2020	01	Addition of EC Rep, CE mark, notified body number
3.	09/02/2021	02	Change of Notified body number, and correction of postal code
4.	01/06/2023	03	IFU translated to Dutch, French & German Languages, No changes in IFU itself.
5.	19/03/2024	04	Additional Symbols added, "Notice" added
6.	19/06/2024	05	IFU translated to Portuguese, & Spanish Languages, No changes in IFU itself.

Sample of the Information Release document:

 <p>alspl</p> <p>Changing Lives For The Better....</p>	<b>Advanced LifeSciences Pvt. Ltd.</b>	<b>Document No.-</b> QSP/SOP/051		
		<b>Page No.</b>	Page 2 of 22	
	<b>Product Information Release</b>		<b>Effective Date</b>	01/06/2023
			<b>Revision No.</b>	01
			<b>Applicable Document:</b>	QSP/SOP/05

**Surgewire™ and Hydrowire™**  
**Diagnostic Guidewires**  
**Voerdraden**  
**Fils Guides**  
**Führungsdrähte**  
**Fio-Guia**  
**Alambre Guía**

**Table of Contents**

<a href="#"><u>ENGLISH (en)</u></a> .....	3
<a href="#"><u>DUTCH (nl)</u></a> .....	5
<a href="#"><u>FRENCH (fr)</u></a> .....	8
<a href="#"><u>GERMAN (de)</u></a> .....	10
<a href="#"><u>PORTUGUESE (pt)</u></a> .....	13
<a href="#"><u>SPANISH (es)</u></a> .....	16

 0297

 alspl <small>Changing Lives For The Better....</small>	<b>Advanced LifeSciences Pvt. Ltd.</b>	<b>Document No.-</b> QSP/SOP/051		
		<b>Page No.</b>	Page 3 of 22	
	<b>Product Information Release</b>		<b>Effective Date</b>	01/06/2023
			<b>Revision No.</b>	01
			<b>Applicable Document:</b>	QSP/SOP/05

**ENGLISH (en)**

## Instructions for Use

### Description

A guidewire is a wire, coil or spring used as a guide for placement of a larger device such as a catheter. The guidewire usually consists of a core wire made of different metallic materials like stainless steel or nitinol, a spring / coil wire which act as a coating, or other forms of coating like PTFE, hydrophilic coatings or other materials based on the user's need or preference. The device comes in varying sizes, lengths and designs (materials and tips).

Advanced LifeSciences Pvt. Ltd.'s (ALSPL) guidewires are packaged in a plastic dispenser which is fitted with a luer hub at one end. This helps in flushing the device before use as recommended in the Instruction for Use below.

### Indications for Use:

ALSPL's guidewires are intended to facilitate percutaneous vessel entry or placement of devices during diagnostic and interventional procedures, typically using the Seldinger technique. The guidewires are used for urology, gastroenterology or peripheral vascular applications

### Contraindications:

- Cerebral and central circulatory vasculature.
- As the guidewire is only an accessory for procedures, the contraindications are limited to the contraindications for the respective procedures

### Warnings and Precautions:

1. Please read all packaging insert warnings, precautions, and instructions before use. Failure to do so may result in severe patient injury or death.
2. Warning: Sterile by EO, Single Use: Do not reuse, reprocess or re-sterilize under any circumstance. Blood and other foreign material can collect in the inner parts of the guidewire which cannot be removed; therefore, guidewires are for single use only.
3. Warning: Do not use if the packaging is damaged or opened unintentionally before use. This will lead to loss of sterility.
4. Precaution: Aseptic techniques must be used wherever possible, especially during removal from packaging.
5. Precaution: The guidewire must be carefully and gently used with fluoroscopic guidance to facilitate the complex movement within the body.
6. Warning: Rx Only. The guidewire must be used by trained and experienced personnel, who are well versed in catheterization techniques, anatomical landmarks, and potential complications.
7. Precaution: The product must be stored in a dry environment, away from sunlight.
8. Precaution: PTFE and Hydrophilic guidewires should not be withdrawn with metallic needles as those can damage or separate the coating of the wires. Metallic needles should be replaced with introducer sheaths, catheters, or vessel dilators as soon as possible.
9. Warning: Guidewires contain metal and should not be used during MRI procedures.
10. Warning: If resistance is felt during use, the manipulation of the guidewire, and the accompanying device (catheter) should be stopped, and the cause of resistance should be

 <p>alspl</p> <p>Changing Lives For The Better...</p>	<b>Advanced LifeSciences Pvt. Ltd.</b>	<b>Document No.-</b> QSP/SOP/051		
		<b>Page No.</b>	Page 4 of 22	
	<b>Product Information Release</b>		<b>Effective Date</b>	01/06/2023
			<b>Revision No.</b>	01
			<b>Applicable Document:</b>	QSP/SOP/05

determined if possible. If removal is required, the guidewire and the accompanying device should be removed together as a unit.

11. Precaution: Guidewires may lead to thrombosis, as is the case with most interventional products. The user is responsible for taking adequate measures to prevent such an occurrence.
12. Warning: Do not use if the product is kinked or damaged in any way.
13. In the event that there is a serious incident that has occurred in relation to the device, please report such an occurrence to the company (contact details at the bottom) immediately, and to the competent authority in your state.
14. This product is only designed, intended for use as indicated.
15. Please do not alter the guidewire in any way during the procedure.

#### **Additional warnings and precautions for Hydrophilic Guidewires**

16. Precaution: If the guidewire is to be reinserted in a patient after use, it should be cleaned in a heparinised saline solution, or wiped with a gauze dipped in heparinised saline solution. Use of alcohol, antiseptic solutions must be avoided as they can damage the hydrophilic coating
17. Precaution: Only use the torquers provided with the Hydrophilic guidewires for manipulation. Use of metallic torquers can damage the hydrophilic coating of the wires.

#### **Instruction for Use:**

##### **Flushing the guidewire before use:**

1. Fill a syringe (preferably with a luer lock) with saline or heparinised saline solution
2. Attach the filled syringe to the luer hub of the guidewire dispenser
3. Inject the solution into the dispenser
4. Detach the syringe from the dispenser
5. Re-flush the guidewire if needed, or remove the guidewire from the dispenser

##### **Guidewire insertion:**

6. The guidewire can be inserted using the well-known Seldinger technique, where the puncture site is punctured using an introducer needle, and the guidewire is introduced through the needle
  - 6.1. When this technique is used, the needle must be removed as soon as possible and replaced with an appropriate dilator or introducer sheath to avoid damage to the guidewire
7. The guidewire can also be inserted in a catheter or other interventional device already in the body
  - 7.1. In such a case, the guidewire lumen of the accompanying device must be flushed before insertion of the guidewire in the lumen
  - 7.2. At least 5cm of the guidewire should always be extended out of the accompanying device
8. ALSPL guidewire dispensers have straighteners at the tip to aid in insertion of 'J' tipped, or angled tipped guidewires

##### **Guidewire Disposal:**

9. Guidewires can be a biological hazard if disposed incorrectly. They should be disposed using recognised medical procedures and industry best practices. Local laws and regulations should be taken into account during disposal

#### **Potential Complications / Adverse Effects**

Guidewires are accessories used in various medical procedures. The potential complications or adverse effects caused during use are limited to the procedure being performed.

#### **Warranty and Limitations**

 Changing Lives For The Better...	<b>Advanced LifeSciences Pvt. Ltd.</b>	<b>Document No.-</b> QSP/SOP/051		
		<b>Page No.</b>	Page 5 of 22	
	<b>Product Information Release</b>		<b>Effective Date</b>	01/06/2023
			<b>Revision No.</b>	01
			<b>Applicable Document:</b>	QSP/SOP/05

Guidewires are sold in 'as is' condition. The entire risk as to the quality and performance of the product is with the buyer. ALSPL disclaims all warranties, expressed or implied, with respect to its products, including but not limited to, any implied warranty of merchantability or fitness for a particular purpose. ALSPL shall not be liable to any person for any medical expenses or any direct or consequential damages resulting from the use of any product or caused by any defect, failure, or malfunction of any product, whether a claim for such damages is based upon warranty, contract, tort, or otherwise. No person has any authority to bind ALSPL to any representation or warranty with respect to its products.

**Notice:** Any serious incident which occurs in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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 Changing Lives For The Better...	<b>Advanced LifeSciences Pvt. Ltd.</b>	<b>Document No.-</b> QSP/SOP/051		
		<b>Page No.</b>	Page 6 of 22	
	<b>Product Information Release</b>		<b>Effective Date</b>	01/06/2023
			<b>Revision No.</b>	01
			<b>Applicable Document:</b>	QSP/SOP/05

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	<b>Advanced LifeSciences Pvt. Ltd.</b>	<b>Document No.-</b> QSP/SOP/051		
		<b>Page No.</b>	Page 7 of 22	
	<b>Product Information Release</b>	<b>Effective Date</b>		01/06/2023
		<b>Revision No.</b>		01
		<b>Applicable Document:</b>		QSP/SOP/05

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 Changing Lives For The Better...	<b>Advanced LifeSciences Pvt. Ltd.</b>	<b>Document No.-</b> QSP/SOP/051		
		<b>Page No.</b>	Page 8 of 22	
	<b>Product Information Release</b>		<b>Effective Date</b>	01/06/2023
			<b>Revision No.</b>	01
			<b>Applicable Document:</b>	QSP/SOP/05

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 Changing Lives For The Better...	<b>Advanced LifeSciences Pvt. Ltd.</b>	<b>Document No.-</b> QSP/SOP/051	
		<b>Page No.</b>	Page 9 of 22
<b>Product Information Release</b>	<b>Effective Date</b>		01/06/2023
	<b>Revision No.</b>		01
	<b>Applicable Document:</b>		QSP/SOP/05

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 <small>Changing Lives For The Better...</small>	<b>Advanced LifeSciences Pvt. Ltd.</b>	<b>Document No.-</b> QSP/SOP/051	
		<b>Page No.</b>	Page 10 of 22
<b>Product Information Release</b>		<b>Effective Date</b>	01/06/2023
		<b>Revision No.</b>	01
		<b>Applicable Document:</b>	QSP/SOP/05

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	<b>Advanced LifeSciences Pvt. Ltd.</b>	<b>Document No.-</b> QSP/SOP/051		
		<b>Page No.</b>	Page 11 of 22	
	<b>Product Information Release</b>		<b>Effective Date</b>	01/06/2023
			<b>Revision No.</b>	01
			<b>Applicable Document:</b>	QSP/SOP/05

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	<b>Advanced LifeSciences Pvt. Ltd.</b>	<b>Document No.-</b> QSP/SOP/051		
		<b>Page No.</b>	Page 12 of 22	
	<b>Product Information Release</b>		<b>Effective Date</b>	01/06/2023
			<b>Revision No.</b>	01
			<b>Applicable Document:</b>	QSP/SOP/05

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 <small>Changing Lives For The Better...</small>	<b>Advanced LifeSciences Pvt. Ltd.</b>	<b>Document No.-</b> QSP/SOP/051	
		<b>Page No.</b>	Page 13 of 22
<b>Product Information Release</b>		<b>Effective Date</b>	01/06/2023
		<b>Revision No.</b>	01
		<b>Applicable Document:</b>	QSP/SOP/05

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 <p>alspl</p> <p>Changing Lives For The Better...</p>	<b>Advanced LifeSciences Pvt. Ltd.</b>	<b>Document No.-</b> QSP/SOP/051		
		<b>Page No.</b>	Page 14 of 22	
	<b>Product Information Release</b>		<b>Effective Date</b>	01/06/2023
			<b>Revision No.</b>	01
			<b>Applicable Document:</b>	QSP/SOP/05

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 <p>alspl</p> <p>Changing Lives For The Better...</p>	<b>Advanced LifeSciences Pvt. Ltd.</b>	<b>Document No.-</b> QSP/SOP/051		
		<b>Page No.</b>	Page 15 of 22	
	<b>Product Information Release</b>		<b>Effective Date</b>	01/06/2023
			<b>Revision No.</b>	01
			<b>Applicable Document:</b>	QSP/SOP/05

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	<b>Advanced LifeSciences Pvt. Ltd.</b>	<b>Document No.-</b> QSP/SOP/051		
		<b>Page No.</b>	Page 16 of 22	
	<b>Product Information Release</b>	<b>Effective Date</b>	01/06/2023	
		<b>Revision No.</b>	01	
		<b>Applicable Document:</b>	QSP/SOP/05	

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 <p>alspl</p> <p>Changing Lives For The Better...</p>	<b>Advanced LifeSciences Pvt. Ltd.</b>	<b>Document No.-</b> QSP/SOP/051		
		<b>Page No.</b>	Page 17 of 22	
	<b>Product Information Release</b>		<b>Effective Date</b>	01/06/2023
			<b>Revision No.</b>	01
			<b>Applicable Document:</b>	QSP/SOP/05

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 <small>Changing Lives For The Better...</small>	<b>Advanced LifeSciences Pvt. Ltd.</b>	<b>Document No.-</b> QSP/SOP/051	
		<b>Page No.</b>	Page 18 of 22
<b>Product Information Release</b>		<b>Effective Date</b>	01/06/2023
		<b>Revision No.</b>	01
		<b>Applicable Document:</b>	QSP/SOP/05

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

 Changing Lives For The Better....	<b>Advanced LifeSciences Pvt. Ltd.</b>	<b>Document No.-</b> QSP/SOP/051		
		<b>Page No.</b>	Page 19 of 22	
	<b>Product Information Release</b>		<b>Effective Date</b>	01/06/2023
			<b>Revision No.</b>	01
			<b>Applicable Document:</b>	QSP/SOP/05

	Fecha de manufactura
	Do not re-use Geen hergebruik ne pas réutiliser Nicht wiederverwenden Não reutilize No reutilizar
	Do not use if package is damaged Geen gebruik als de verpakking is beschadigd ne pas utiliser si l'emballage est endommagé Nicht verwenden, wenn die Verpackung beschädigt ist Não use se a embalagem estiver danificada No utilizar si el paquete está dañado.
	Sterilized using ethylene oxide Gesteriliseerd met ethyleenoxide Stérilisé à l'oxyde d'éthylène Mit Ethylenoxid sterilisiert Esterilizado com óxido de etileno Esterilizado con óxido de etileno.
	Keep away from sunlight Geen blootstelling aan zonlicht Tenir à l'écart de la lumière du soleil Von Sonnenlicht fernhalten Mantenha longe da luz solar Mantener alejado de la luz solar
	Batch code Batch code Code du lot Chargencode Código de lote Código de lote
	Catalogue Number Catalogus nummer Numéro de catalogue Katalog-Nummer Catálogo de número Número de catalogo
	Manufacturer Producent Fabricant Hersteller Fabricante Fabricante
	Use-by date Vervaldatum Date de péremption Verfallsdatum Data de validade Utilizar por fecha

 Changing Lives For The Better...	<b>Advanced LifeSciences Pvt. Ltd.</b>	<b>Document No.-</b> QSP/SOP/051		
		<b>Page No.</b>	Page 20 of 22	
	<b>Product Information Release</b>		<b>Effective Date</b>	01/06/2023
			<b>Revision No.</b>	01
			<b>Applicable Document:</b>	QSP/SOP/05

	Keep dry Hou droog garder au sec Trocken aufbewahren Manter seco Mantener seco
	Temperature Limit(5-45°C) Toegelaten temperatuur(5-45°C) limite de temperature(5-45°C) Temperatur-Grenzwert(5-45°C) Limite de temperatura (5-45°C) Temperatura límite(5-45°C)
	Authorised representative in the European union Gemachtigd vertegenwoordiger in de Europese Unie Représentant autorisé dans l'Union européenne Bevollmächtigter Vertreter in der Europäischen Union Representante autorizado na União Europeia Representante autorizado en la unión europea
	Non-Pyrogenic Niet-pyrogeen Non-pyrogène Nicht pyrogene Não Pirogênico No pirogénico
	Prescription use only Alleen gebruik op recept Utilisation sur ordonnance uniquement verschreibungspflichtig verwenden Nur Somente uso com receita Sólo para prescripción
	Medical device Medisch apparaat Dispositif médical Medizinisches Gerät Aparelho médico Dispositivo médico
	Do not re-sterilize Niet opnieuw steriliseren Ne pas restériliser Nicht erneut sterilisieren Não reesterilize No reesterilizar
	Single Sterile Barrier System Enkelvoudig steriel barrièresysteem Célibataire Système de barrière stérile Einzelnes Sterilbarrieresystem Sistema de barreira estéril solteira Sistema de barrera estéril única

 Changing Lives For The Better...	<b>Advanced LifeSciences Pvt. Ltd.</b>	<b>Document No.-</b> QSP/SOP/051		
		<b>Page No.</b>	Page 21 of 22	
	<b>Product Information Release</b>		<b>Effective Date</b>	01/06/2023
			<b>Revision No.</b>	01
			<b>Applicable Document:</b>	QSP/SOP/05

	<p>Single Sterile barrier system with protective packing outside          Enkel steriel barrièresysteem met beschermende verpakking aan de buitenkant          Système de barrière stérile Célibataire avec emballage de protection à l'extérieur          Einzel-Sterilbarrieresystem mit außenliegender Schutzverpackung          Sistema de barreira estéril único com embalagem protetora externa          Sistema de barrera estéril soltera con embalaje protector exterior</p>
	<p>Double Sterile Barrier System          Dubbel steriel barrièresysteem          Système de double barrière stérile          Doppeltes Sterilbarrieresystem          Sistema de barreira dupla estéril          Sistema de doble barrera estéril</p>
	<p>Single sterile barrier system with protective packing inside          Enkelvoudig steriel barrièresysteem met beschermende verpakking aan de binnenkant          Système de barrière stérile Célibataire avec emballage de protection à l'intérieur          Einzelnes Sterilbarrieresystem mit Schutzverpackung im Inneren          Sistema de barreira estéril solteira com embalagem protetora interna          Sistema de barrera estéril soltera con embalaje protector interior</p>

 <b>Manufactured By:</b> Advanced LifeSciences Pvt. Ltd. D-22 Okhla Industrial Area Phase – 1, New Delhi – 110020, India. Website : <a href="http://www.alspl.com">www.alspl.com</a> , E-mail: <a href="mailto:customerservice@alspl.com">customerservice@alspl.com</a> Phone: +91 9818237529	<b>EC REP</b>
	Name: CMC Medical Devices & Drugs S.L. Address: C/ Horacio Lengo N° 18, CP 29006, Málaga, Spain
	<b>UK REP</b>
	Name: BVM Medical Address: BVM House Trinity Lane Hinckley, Leicestershire LE10 OBL, United Kingdom.

IFU/14/01

Revision: 05 19/06/2024

Regulatory Review:

Sr. No.	Scope	Check	Applicable Section(s)	Comments
1.	Indian MDR 2017	<input checked="" type="checkbox"/>	-	-
2.	EU MDR 2017/745	<input checked="" type="checkbox"/>	Annex 1, Section 23.4	-
3.	US FDA	<input type="checkbox"/>	-	-
4.	IS / EN ISO 13485	<input checked="" type="checkbox"/>	-	-
5.	Any Other, please mention	<input type="checkbox"/>	-	-

 Changing Lives For The Better...	<b>Advanced LifeSciences Pvt. Ltd.</b>		<b>Document No.-</b> QSP/SOP/051	
			<b>Page No.</b>	Page 22 of 22
	<b>Product Information Release</b>		<b>Effective Date</b>	01/06/2023
			<b>Revision No.</b>	01
			<b>Applicable Document:</b>	QSP/SOP/05

<b>Release product info</b>	<b>Responsible</b>	<b>Date</b>	<b>Name/ Function (print)</b>	<b>Signature</b>
Created	Creator	19/06/2024	Rishabh Pandey	
CC specifications checked	Production Department	19/06/2024	Ajit Kumar Sahu	
Content / specifications checked	QA/QC Department	19/06/2024	Kusum Gupta	
Release: Specifications and regulatory requirements reviewed	RA	19/06/2024	Shubhankar Aggarwal	